

Product Information

Other Animal Sera

Product Description

Other animal sera are used in various applications, e.g. autologous cell culture systems or immune-histological methods.

Porcine serum is often used in cultures that proliferate mycoplasma or viruses. Goat and Sheep Serum are used in immunohistochemistry and to eliminate non-specific binding reactions.

Applications

- Autologous cell culture systems
- Immunohistochemistry
- Immunoblotting
- Immunofluorescence
- Immunoprecipitation

Test Parameter

- pH
- Osmolality
- Endotoxin
- Total protein
- Hemoglobin
- Mycoplasma
- Sterility

Storage: Store at $\leq -15^{\circ}\text{C}$

Product	Volume	Cat. No.
Goat Serum	100 ml	CG100-1010
	500 ml	CG500-1010
Porcine Serum	100 ml	CP100-1010
	500 ml	CP500-1010
Sheep Serum	100 ml	CS100-1010
	500 ml	CS500-1010
Mouse Serum	100 ml	CM100-1010
Rat Serum	100 ml	CRT100-1010
Rabbit Serum	100 ml	CRB100-1010

Precaution and Disclaimer

This product is for research use only. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.