

Product Information

Medium 199 with Earle's Salts, with L-Glutamine

Cat. No.: E0500-260 Volume: 500 ml

Product Specifications

CO ₂ concentration, optimum	5.0 %
pH	7.0 – 7.5
Osmolality	280 – 340 mOsm/kg
Colour and Appearance	Clear red orange solution
Endotoxin	< 1 EU/ml (< 0.1 ng/ml)
Cell Culture	Tested
Sterility	Tested
Storage	+2°C to +8°C protected from light

Formulation

Components	Concentration mg/l	Molecular Weight g/mol	Concentration mM
Amino Acids			
L-Alanine	50.00	89.09	0.56
L-Arginine HCl	70.00	210.66	0.33
L-Aspartic Acid	60.00	133.10	0.45
L-Cysteine HCl H ₂ O	0.11	176.00	0.0006
L-Cystine 2HCl	26.00	240.30	0.11
L-Glutamic Acid	133.60	147.13	0.91
L-Glutamine	100.00	146.15	0.68
Glycine	50.00	75.07	0.67
L-Histidine HCl H ₂ O	21.88	209.63	0.10
L-Hydroxy-L-Proline	10.00	131.00	0.08
L-Isoleucine	40.00	131.17	0.30
L-Leucine	120.00	131.17	0.91
L-Lysine HCl	70.00	182.65	0.38
L-Methionine	30.00	149.21	0.20
L-Phenylalanine	50.00	165.19	0.30
L-Proline	40.00	115.13	0.35
L-Serine	50.00	105.09	0.48
L-Threonine	60.00	119.12	0.50
L-Tryptophan	20.00	204.23	0.10
L-Tyrosine 2Na 2H ₂ O	57.66	225.15	0.26
L-Valine	50.00	117.15	0.43

Components	Concentration mg/l	Molecular Weight g/mol	Concentration mM
------------	--------------------	------------------------	------------------

Product Information

Medium 199 with Earle's Salts, with L-Glutamine

Cat. No.: E0500-260 Volume: 500 ml

Vitamins				
Ascorbic Acid	0.056	176.12	3.18E-4	
D-Biotin	0.01	244.31	4.09E-5	
Choline chloride	0.50	139.62	0.0036	
D-Calcium Pantothenate	0.01	476.53	2.10E-5	
Folic Acid	0.01	441.40	2.27E-5	
Menadione Sodium Bisulfite	0.016	276.24	5.79E-5	
myo-Inositol	0.05	180.16	2.78E-4	
Nicotinamide	0.025	122.12	2.05E-4	
Nicotinic Acid	0.025	123.11	2.03E-4	
Paraaminobenzoic Acid	0.05	137.14	3.65E-4	
Pyridoxal HCl	0.025	203.62	1.23E-4	
Pyridoxine HCl	0.025	205.64	1.22E-4	
Riboflavin	0.01	376.36	2.66E-5	
Thiamine HCl	0.01	337.27	2.96E-5	
Vitamin A (Acetate)	0.14	328.49	4.26E-4	
Vitamin D2 (Ergocalciferol)	0.10	396.65	2.52E-4	
(±)-α-Tocopherol Phosphate 2Na	0.01	554.65	1.80E-5	
Inorganic Salts				
Calcium Chloride 2H ₂ O	265.00	147.01	1.80	
Ferric Nitrate 9 H ₂ O	0.72	404.00	0.0018	
Magnesium Sulfate	97.67	120.37	0.81	
Potassium Chloride	400.00	74.55	5.37	
Sodium Bicarbonate	1250.00	84.01	17.88	
Sodium Chloride	6800.00	58.44	116.36	
Sodium Phosphate Monobasic Anhydrous	122.00	119.98	1.02	
Other Components				
2-Deoxy-D-Ribose	0.50	134.13	0.004	
Adenine Sulfate	10.00	184.17	0.05	
Adenosine 5'-Phosphate	0.24	365.24	6.57E-4	
Adenosine 5'-Triphosphate 2Na	1.00	551.14	0.002	
Cholesterol	0.20	386.65	5.17E-4	
D-Glucose	1000.00	180.16	5.55	
L-Glutathione (reduced)	0.05	307.32	1.63E-4	
Guanine	0.30	151.13	0.002	
Hypoxanthine	0.30	136.11	0.002	
Components	Concentration mg/l	Molecular Weight g/mol	Concentration mM	
Phenol Red Sodium Salt	21.3	376.36	0.057	

Product Information

Medium 199 with Earle's Salts, with L-Glutamine

Cat. No.: E0500-260 Volume: 500 ml

D-Ribose	0.50	150.13	0.003
Sodium Acetate	50.00	82.00	0.61
Thymine	0.30	126.11	0.002
Tween® 80	20.00	-	Infinity
Uracil	0.30	112.09	0.003
Xanthine Na	0.34	174.09	0.002

Tween is a registered trademark of Croda International PLC.

Precautions and Disclaimer

This product is for research use only. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.